

CURIE

Versatility that Delivers



Enlightening your creations

CURIE is equipped with a sophisticated UV-light source management system that enables users to do more in less time and space. Ideal for dental professionals who take bio-compatibility seriously and care about patient safety



Efficient Curing Space

Perfect curing chamber space for Dental Clinics and Labs with rigorous workflows.



Future-proof

Equipped with an SD card reader for updating material settings and post-curing times.



Accuracy and consistency

Preloaded optimal curing times database available at the turn of a dial and touch of a button.



Premium Performance

Smart UV-light source management system for maximized energy efficiency and the ideal light strength to prevent any deformation.

TECHNICAL SPECIFICATIONS



Power Input

100 - 240 VAC / 50-60 Hz



LED Power Output

70 W



Wavelength of LED

365, 385, 405 nm



Curing Time

1 sec - 20 min



Dimensions

L 20 x W 22 x H 16 cm (Max)
L 7.87 x W 8.66 x H 6.30 inch (Max)



Weight

3 Kg
6.61 lbs



Chamber Dimensions

L 14.2 x W 14.5 x H 8 cm (Max)
L 7.87 x W 8.66 x H 6.30 inch (Max)

PACKAGE INCLUDES

CURIE machine	1 Unit
Glass plate (Inside the machine)	1 Pcs
Power cord	1 Pcs
SD card	1 Pcs
Card reader	1 Pcs
Quick start guide	1 Pcs